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

- | | |
|-------------------------------------|-----------------------------------|
| 1. Strap for transport | 5. Clamp |
| 2. Elastomeric pump | 6. Filter |
| 3. Check-Valve for filling with cap | 7. Luer lock male connector |
| 4. Extension line | 8. Luer lock male cap with filter |

PRODUCT FEATURES AND USE:

Single use medical device for the clinical treatment of pain. The elasticity of the silicone allows a constant control of flow through the catheter capillary ensuring a safe and effective infusion during therapy.

PURPOSE:

The product can be used in continuous infusion for the clinical treatment of pain, for cancer therapy, for infusion of antibiotics and in the treatment of thalassemias therapies.

PRECAUTIONS AND CONTRAINDICATIONS:

- Sterile, non-toxic, non-pyrogenic device. The expiry date is 5 years from the date of sterilization and refers to the product distributed and stored as indicated below.
 - Use only if the package is not damaged or tampered with. Do not use if the caps are not included or are disconnected.
 - Store in a clean place at room temperature, away from heat and moisture.
 - Device for single use. Do not fill twice. Do not reuse. In case of re-use it is possible a loss of functionality of the system and risk of infection to the patient.
 - Discard after use according to local regulations.
 - The elastomeric pump is not intended for the administration of blood, haemoderivates or TPN lipid emulsions.
 - Respect the nominal volume of filling indicated on the label of the device. The filling volume is given by the nominal volume plus residual volume (no more than 10%). Do not fill the tank with a volume less than or greater than the nominal one.
 - Different filling volumes different from the nominal ones can lead to significant fluctuations in the delivered flow.
 - The flow may vary according to the viscosity and concentration of the drug infusion and temperature of the environment and the patient. In particular, the models with PCA should be worn outside of the clothing and stored at room temperature.
 - The speed of the flow is optimal when the elastomeric reservoir and the distal luer lock are placed at the same height.
 - Infusion time may increase significantly after a prolonged storage period: It is recommended to enable the device temperature to stabilize before use.
 - Before introducing the drug, insert the saline solution to dilute the drug itself: this precaution avoids the administration of an undiluted drug.
 - Administer the prescribed volume of medication: overdose could alter the outcome of the therapy.
 - When using a catheter, refer to the instructions for use provided by the manufacturer of the catheter. The length, diameter and position of the catheter can affect the flow. Do not use a catheter with a diameter less than 22 gauge (3French).
- Free of latex components and phthalates (DEHP). The models in PVC may be incompatible with some pharmacological solutions: refer to instruction for use of the drug for relevant information. The device is suitable for home therapy.
- Anti-UV models are suitable for photosensitive solutions. The device is not suitable for neonatology.

INSTRUCTIONS FOR USE:

The elastomeric pump is filled using a syringe with luer. It is recommended to use a filter during filling.

- Check the clarity of the solution before administration, do not use if it is not clear.
 - Check the integrity of the package and take the elastomeric pump just before use.
 - Use an aseptic technique throughout the procedure.
 - Remove the cap from the filling valve, close the line with the clamp,
 - Secure the syringe to the filling valve without overtightening and introduce the drug keeping the system in vertical position.
 - The filling valve for the addition of the solution is a check-valve: if it is necessary to reduce the solution already introduced in the elastomer, take the special individually sterile packaged adapter, connect the adapter to the syringe, insert in the check-valve to activate the drug withdrawal.
 - During the infusion of the drug, make sure there is no air in the syringe. If there is air in the elastomeric silicone, rotate the syringe and suck up all the air (if the air is included, the elastomeric silicone automatically releases the air in maximum two hours)
 - After introducing the desired volume, remove the syringe and replace the protective cap on the filling check-valve. Open the clip.
- The beginning of dripping from the end part is an indication of the correct functioning of the device
- In the models with bolus PCA, press the button once or twice to activate the infusion of the liquid in the line. After checking the presence of fluid at the end of the line, close the clip and the male luer lock fitting at the end of the line to use the device next time.
- The PCA button is the function button that the patient can use to control the additional medicine when he is in a continuous infusion position. By pressing the PCA button, the patient can add a limited amount of medicine according to the medical disposition, as needed.
 - If the medicine does not flow, connect a luer adapter or a 3-way stopcock to the distal end of the elastomer. Connect a syringe to the other side and keep it on the suction until the air is completely removed.
 - Remove the syringe, check that the liquid begins to flow out of the fitting and replace the protective cap on the distal part of the system.
 - Label the elastomer with preparation date, patient name, drugs used and their dosage.
 - Check the clarity of the solution both before positioning the system and in the following days.
 - Mark the infusion start level with a mark of indelible marker to have a clear reference for the functioning checks.
 - Connect the system to a suitable device for administering to the patient.
 - If, during administration, leaks occur at any point on the line, stop the infusion cycle, check the source of leak and replace the entire system.

The flow can be slightly faster during the first 1-2 hours of use. This is due to the physical characteristics of the silicone elastomer. The average flow rate per nominal flow value is $\pm 15\%$, determined with distilled water, at room temperature (23 ± 2) °C, relative humidity ($50 \pm 5\%$) and atmospheric pressure 86 kpa ~ 106KPa.

The flow may vary due to:

1) Volume of filling:

When used, add the liquid to the nominal volume, a too high or too low filling volume will cause an imprecise flow of infusion.

2) Viscosity and / or concentration of the drug:

The flow rate has been determined with purified water or distilled water, the infusion of an excessively viscous liquid will cause a slower flow.

3) Temperature:

The flow rate has been determined at the temperature of ($23 + 2$) °C, the flow velocity will be rapid when the temperature of use is higher, in contrast at lower temperatures, the flow will be slower.

4) Atmospheric pressure:

The product flow rate has been determined under standard atmospheric pressure conditions, the flow rate of the product will be rapid when the use temperature is higher, vice versa will be slower.

5) Infusion level:

The normal use of the product is for infusion, the suspension of the infusion can lead to a faster flow.

6) Conservation:

The product should be used promptly after filling, the flow delivered will be slower after long storage.

VARIATION OF INFUSION SPEED (only for models with Multirate system)

FOR ITEMS WITH FLOW REGULATION SYSTEM M1: Follow the directions described above and:

- Insert the flow regulator key, in plastic poche, inside the "control module" and change the speed until the one you want to start filling the elastomer.

- make the priming of the line.

- remove the key and store it. **WARNING:** Be sure not to leave the key between two flow values settings: the system may give errors.

FOR ITEMS WITH FLOW REGULATION SYSTEM M2 AND M3: Follow the directions described above and:

- rotate the flow regulator key to the maximum value allowed

- make the priming of the line.

- rotate the flow regulator key to the desired value, then remove and store it. **WARNING:** Be sure not to leave the key between two flow values settings: the system may give errors.

ACTIVATION OF BOLUS (only for dedicated models)

The button for the bolus dose for PCA allows the patient to independently manage additional infusion, 0.5 ml of drug during continuous infusion. If necessary, press the button for the bolus to PCA. The following bolus is available after 15 minutes.

Available high-flow models

For all models the amount of solution actually administered is > 90%

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