Multimedical s.r.l.
Via Guido Rossa, 71 - 46019 Viedana (MN)
Tel. 0039-0375-765882 Fax 0039-0375-765885
www.multimedical.it

CE 0123

MD

STERILE ISO

10402-2020/07

I - AGO DI HUBER

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10402-2020/07

GB-HUBER NEEDLE

INFUSION SYSTEM FOR INSERTION INTO THE SEPTUM OF SUBCUTANEOUS IMPLANTED PORT

90° needle. The needle size and the length of the tube connection are shown on the product label. ATOMIC AND NON PYROGENIC. THE PRODUCT DOES NOT CONTAIN LATEX AND PHTHALATES.

WARNINGS:
Use only if package is not damaged and tampered with. Do not use the device if the caps are missing or detached. Keep in clean and dry place, away from moisture, direct sunlight and heat source. The device is for single use only. Do not re-use. Reuse can use a loss of device functionality and a possible risk of patient infection. Dispose of the product after its use according to the current regulations. Do not use with drugs and/or substances that are incompatible with PVC. The device is not suitable for use in neonatology. Not suitable for photosensitive solutions. The connection fittings of this device comply with ISO 80369-7. The device must be thrown into a biohazard container. Dispose according to national / international standards and hospital protocols. Use is only allowed to qualified personnel.

Notice: In the event of an incident that has occurred in relation to the device, report the event to the manufacturer and the Competent Authority.

INSTRUCTIONS FOR USE:

- Open the package immediately before use.
- Practice a local disinfection of the implantation area, wash and disinfect hands;
- connect the infusion disposal to a syringe and expel the air from the set;
- disinfect the area of puncture.
- Localize the center of the port by palpation maintaining firmly the port between 2 fingers;
- insert needle perpendicular by holding the special fins and advance the needle to the base of the port.
- confirm correct needle placement in the port reservoir by aspiration of blood before infusion of any substances.
- Fix the needle on the patient with transparent medication and start with planned therapy.
- If during the administration should occur leak, replace the entire system.
- In the case of impromptu infusions of medications, use the injection point of the Y-fitting, where present (do not use needles of diameter greater than 23G)
- Only use the correct fit size of non-coring Huber needle. If needle is too short, the reservoir will not be entered, and it will not be possible to aspirate and administer medications. If the needle is too long, it will protrude and potentially dislodge. Be careful to insert the needle completely: an incomplete insertion can cause partial occlusion of the needle by the septum of the port with the risk of difficult infusion and / or dislocation of the needle. Do not tilt, rock, rotate or pull on needle when accessing port.
- Potential complications:
- Infection of the site of needle puncturing,
- Needle dislocation and drug extravasation,
- Air embolism.

Drug	Viscosity	Diameter
NaCl 0,9%	Low	G24-G22
Plasma	Low	G22
Lipidic emulsion 20%	Low	G22-G20
Lipidic emulsion 30%	Medium	G22-G20
Cytotoxic		
Glucose 30%	High	G20-G19
Glucose 50%	High	G19
Haemo-derivates		
Blood samples		

Needle length	Directions for use
25 mm	Standard size port obese patients
20 mm	Standard port and low profile Normal size patients
15 mm	Low profile Port Thin patients
12 mm	Baby port or brachial port surface facilities children

Multimedical s.r.l.
Via Guido Rossa, 71 - 46019 Viedana (MN)
Tel. 0039-0375-765882 Fax 0039-0375-765885
www.multimedical.it

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